



# Manual

## Pre-Installation Guide



Translation

**ELTRA**<sup>®</sup>  
ELEMENTAL ANALYZERS

**Copyright**

© Copyright by  
Eltra GmbH  
Retsch-Allee 1-5  
42781 Haan  
Germany

## Table of Contents

1	Pre-installation guide ELEMENTRAC CS-i .....	4
2	Pre-installation guide ELEMENTRAC CS-d .....	5
3	Pre-installation guide ELEMENTRAC ONH-p, ON-p and OH-p .....	7
4	Pre-installation guide CS 580 .....	9
5	Pre-installation guide CHS 580 .....	10
6	Pre-installation guide CS 580A (Helios) .....	11
7	Pre-installation guide CHS 580A (CHS-Helios) .....	12
8	Pre-installation guide CW 800 .....	13
9	Pre-installation guide CW Multiphase .....	14
10	Pre-installation guide SC 800 .....	15
11	Pre-installation guide H 500 .....	16
12	Pre-installation guide TGA Thermostep, TGA Thermostep ML .....	17
13	Accessories .....	19
13.1	Pre-installation guide Autoloader ELEMENTRAC CS-i/CS-d .....	19
13.2	Pre-installation guide Loader CS 580A (Helios) .....	19
13.3	Pre-installation guide TIC Module .....	19
13.4	Pre-installation guide Gas Purification Furnace GRO18 .....	20
13.5	Pre-installation guide HTF 540 .....	20

## 1 Pre-installation guide ELEMENTRAC CS-i

Following requirements apply, when installing the analyzer:

Carrier gas	Oxygen 99.5% pure; 2 - 4bar
Compressed air	4 – 6bar, oil and grease free
Mains power supply	230VAC ±10%, 50/60Hz; 16A fuse



Analyzer dimension 520 x 840 x 750mm (WHD)

Analyzer weight approx. 150kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 1: Carrier gas tube

### Connections for compressed air:

The tubes supplied together with the analyzer, carry a connector with G¼" inner diameter.




Fig. 2: Compressed air tube

### **NOTICE**

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 2 Pre-installation guide ELEMENTRAC CS-d

Following requirements apply, when installing the analyzer:

Carrier gas                      Oxygen 99.5% pure; 2 - 4bar

Compressed air                4 – 6bar, oil and grease free

**Mains power supply:**

Analyzer:                        230VAC ±10%, 50/60Hz; 16A fuse  
(CEE-Plug 230V, 16A)



Resistance Furnace:        230 VAC ±10%, 50/60 Hz; 20A fuse  
CEE-Plug 230V, 32A



Analyzer dimension        560 x 780 x 600 mm (WHD)

Furnace dimension        340 x 525 x 600 mm (WHD)

Analyzer weight            approx. 140 kg.

Furnace weight            approx. 38 kg.

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.

**Connections for carrier gas:**

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 3: Carrier gas tube

**Connections for compressed air:**

The tubes supplied together with the analyzer, carry a connector with G $\frac{1}{4}$ " inner diameter.




Fig. 4: Compressed air tube

**NOTICE**

**"Internet connection:**

For fast online support connect the analyzer PC to the internet."

### 3 Pre-installation guide ELEMENTRAC ONH-p, ON-p and OH-p

The following prerequisites must be prepared by the customer before installation on site.

Compressed air: 4 – 6bar, oil and grease free

Carrier gas: Pure helium 99.995 %; 2–4 bar, for ONH-p and ON-p  
Pure nitrogen 99.995 %; 2–4 bar, for ONH-p and OH-p  
Pure argon 99.995 %; 2–4 bar; optional

**3-phase version mains voltage supply:**

400 V AC  $\pm$  10 %, L1=32A, L2=30A, L3=0A, 50/60 Hz  
(3L+N+PE, 6h, 32A acc. to IEC 60309)



**1-phase version mains voltage supply:**

230 V AC  $\pm$  10 %; L=63A, 50/60 Hz  
(L+N+PE, 6h, 63A acc. to IEC 60309)



Dimensions of the analyzer: 570 x 770 x 630mm (WHD)

Weight of the analyzer: 161kg

Sound pressure level: 55dB(A)

- The scales need a vibration-free installation site.
- The analyzer needs a sufficiently strong and stable desk as installation site.
- The analyzer requires a connection to an exhausting system or to the outside.

**Connections for carrier gas:**

The supplied tubes carry a connector with G $\frac{1}{4}$ " inner diameter ".




Fig. 5: Carrier gas tube

**Connections for compressed air:**

The tubes supplied together with the analyzer, carry a connector with G $\frac{1}{4}$ " inner diameter.




Fig. 6: Compressed air tube

**Water connection:**

Use tap water for cooling, with 4bar (60psi) pressure,

The tubes supplied together with the analyzer, carry a connector with G $\frac{3}{4}$ " inner diameter.



Fig. 7: Water connection

**NOTICE**

**"Internet connection:**

For fast online support connect the analyzer PC to the internet."



## 4 Pre-installation guide CS 580

Following requirements apply, when installing the analyzer:

Carrier gas	Oxygen 99.5% pure; 2 - 4bar
Mains power supply	230VAC $\pm$ 10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension 560 x 780 x 600mm (WHD)

Analyzer weight approx. 75kg

TIC-module dimension 340 x 525 x 600mm (WHD)

TIC-module weight approx. 28kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G $\frac{1}{4}$ " inner diameter ".




Fig. 8: Carrier gas tube

### NOTICE

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 5 Pre-installation guide CHS 580

Following requirements apply, when installing the analyzer:

Carrier gas	Oxygen 99.5% pure; 2 - 4bar
Mains power supply	230VAC ±10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension 560 x 780 x 600mm (WHD)

Analyzer weight approx. 75kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 9: Carrier gas tube

### NOTICE

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 6 Pre-installation guide CS 580A (Helios)

Following requirements apply, when installing the analyzer:

Compressed air	4 – 6bar, oil and grease free
Carrier gas	Oxygen 99.95 % pure; 2 - 4bar
Mains power supply	230VAC ±10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension	550 x 1000 x 600mm (WHD)
Loader dimension	850 x 350 x 450mm (depending on version) (WHD)
Analyzer weight	approx. 90kg
Loader weight	approx. 55kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 10: Carrier gas tube

### Connections for compressed air:

The tubes supplied together with the analyzer, carry a connector with G¼" inner diameter.




Fig. 11: Compressed air tube

### **NOTICE**

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 7 Pre-installation guide CHS 580A (CHS-Helios)

Following requirements apply, when installing the analyzer:

Compressed air	4 – 6bar, oil and grease free
Carrier gas	Oxygen 99.95 % pure; 2 - 4bar
Mains power supply	230VAC ±10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension	550 x 1000 x 600mm (WHD)
Loader dimension	850 x 350 x 450mm (depending on version) (WHD)
Analyzer weight	approx. 90kg
Loader weight	approx. 55kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 12: Carrier gas tube

### Connections for compressed air:

The tubes supplied together with the analyzer, carry a connector with G¼" inner diameter.




Fig. 13: Compressed air tube

### **NOTICE**

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 8 Pre-installation guide CW 800

Following requirements apply, when installing the analyzer:

Carrier gas	Nitrogen 99.995% pure; 2 - 4bar
Mains power supply	230VAC ±10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension 560 x 780 x 600mm (WHD)

Analyzer weight approx. 75kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 14: Carrier gas tube

### NOTICE

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 9 Pre-installation guide CW Multiphase

Following requirements apply, when installing the analyzer:

Carrier gas	Oxygen / Alternative gas 99.99% pure; 2 - 4bar
Compressed air	4 – 6bar, oil and grease free
Mains power supply	230VAC ±10%, 50/60Hz; 20A fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension 560 x 780 x 600mm (WHD)

Analyzer weight approx. 65 kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 15: Carrier gas tube

### Connections for compressed air:

The tubes supplied together with the analyzer, carry a connector with G¼" inner diameter.




Fig. 16: Compressed air tube

### **NOTICE**

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."

## 10 Pre-installation guide SC 800

Following requirements apply, when installing the analyzer:

Carrier gas: Oxygen 99.995% pure; 2 - 4bar  
Main power supply: 230VAC ±10%, 50/60Hz; 16A fuse.  
(L+N+PE, 6h, 32A acc. to IEC 60309)



Analyzer dimension: 560 x 780 x 600mm (WHD)

Analyzer weight: approx. 65kg

- It is important to install the instrument on a stable place
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G $\frac{1}{4}$ " inner diameter ".




Fig. 17: Carrier gas tube

### NOTICE

If a gas calibration is performed, the above-mentioned carrier gas hose is used.

### NOTICE

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."



## 11 Pre-installation guide H 500

Following requirements apply, when installing the analyzer:

Carrier gas	Nitrogen 99.995% pure; 2 - 4bar
Mains power supply	230VAC ±10%, 50/60Hz; 16A fuse CEE-Plug 230V, 16A



Analyzer dimension 750 x 510 x 600mm (WHD)

Analyzer weight approx. 40kg

- It is important to install the instrument on a stable place.
- The balance should rest on a vibration free support.
- The analyzer requires a connection to an exhausting system or to the outside.

### Connections for carrier gas:

The supplied tubes carry a connector with G¼" inner diameter ".




Fig. 18: Carrier gas tube

### NOTICE

#### "Internet connection:

For fast online support connect the analyzer PC to the internet."



## 12 Pre-installation guide TGA ThermoStep, TGA ThermoStep ML

Following requirements apply, when installing the analyzer:

Carrier gas	Oxygen 99.5% pure; 2 - 4bar Nitrogen 99.5% pure; 2 - 4bar
Compressed air	5 – 6bar, oil and grease free
Mains power supply	230VAC ±10%, 50/60Hz; 25A Fuse (L+N+PE, 6h, 32A acc. to IEC 60309)



### Technical specifications

Exhausting blower throughput – 3,8m<sup>3</sup>/min (maximum)  
 Exhausting air temperature < 70°C (depending on analysis)  
 Electrical power consumption ~ 5,5KW (about 24A / 230VAC)  
 Carousel material – Steel 2.4858  
 Analyzer dimension: 56 x 72 x 53 (WHD closed)  
 Analyzer dimension: 56 x 85 x 92 (WHD opened)  
 Analyzer weight: 65kg

– It is important to install the instrument on a stable place.

### Ventilation

ThermoStep must be vented to an external exhaust.

- The exhausting blower is included in delivery, as well the alumina hose (diameter 100mm)
- To be connected to exhausting system or to outside

### **⚠ CAUTION**

- For connections (gas, electrical and exhausting) the table needs a minimum depth of 90 cm.
- Additional space is needed for gas and exhaust installations in Laboratory.

### Connections for carrier gas for TGA:

The tubes supplied together with the analyzer, carry a connector with G<sup>1</sup>/<sub>4</sub>" inner diameter.




Fig. 19: Carrier gas TGA

**Connections for compressed air:**

The tubes supplied together with the analyzer, carry a connector with G¼" inner diameter.




Fig. 20: Compressed air tube

**NOTICE**

**"Internet connection:**

For fast online support connect the analyzer PC to the internet."

## 13 Accessories

### 13.1 Pre-installation guide Autoloader ELEMENTRAC CS-i/CS-d

Mains power supply      230VAC  $\pm$ 10%, 50/60Hz; 16A fuse  
CEE-Plug 230V, 16A



Loader dimensions:  
Version 36p: 500x250x460mm (WHD)  
Version 130p: 900x500x460 (WHD)

### 13.2 Pre-installation guide Loader CS 580A (Helios)

Loader dimensions:  
Version 36p: 500x250x460mm (WHD)  
Version 130p: 900x500x460 (WHD)

### 13.3 Pre-installation guide TIC Module

Mains power supply      230VAC  $\pm$ 10%, 50/60Hz; 16A fuse  
CEE-Plug 230V, 16A



Dimensions:                      340 x 525 x 600mm (WHD)  
Weight:                            approx. 28kg

### 13.4 Pre-installation guide Gas Purification Furnace GRO18

Mains power supply      230VAC  $\pm$ 10%, 50/60Hz; 16A fuse  
CEE-Plug 230V, 16A



Dimensions:              230 x 300 x 780 mm (WHD)  
Weight:                    approx. 19kg

### 13.5 Pre-installation guide HTF 540

General specifications:  
Furnace temperature:    up to 1550°C  
Mains power supply:    230VAC  $\pm$ 10%, 50/60Hz, 20A fuse  
                                  (L+N+PE, 6h, 32A acc. to IEC 60309)



Dimensions:              340 x 525 x 600mm (WHD)  
Weight:                    approx. 38kg  
–      It is important to install the instrument on a stable place.

---

**ELTRA**<sup>®</sup>  
ELEMENTAL ANALYZERS

**Copyright**

© Copyright by  
Eltra GmbH  
Retsch-Allee 1-5  
42781 Haan  
Germany