

Last Updated 6-11-09

HORIBA PART NUMBER	ALPHA PART NUMBER	DESCRIPTION
900.000.030.001		Teflon Tube dia. 4x6
900.000.040.001		Teflon Tube dia. 6x8
905.110.140.001AN	<u>AR027</u>	Accelerator Particle-type Tungsten 2.5 kg
905.110.140.001AN	<u>AR027</u>	Accelerator Particle-type Tungsten 2 kg
905.110.300.001		Accelerator Pure Iron Chip .5 kg
905.110.300.001AN	<u>AR077</u>	Accelerator Iron 1 kg
905.120.320.001	<u>AR71031</u>	Wire brush, flat type for dust removal
905.130.200.001		Horiba Crucible, without Shelf
905.200.010.001	<u>AR200010</u>	Graphite Crucible box of 1000
905.200.020.001	<u>AR200020</u>	Graphite Crucible, Horiba Style, box of 1000
905.200.080.001	<u>AR241</u>	Silicone Grease for high vacuum
905.200.010.000	<u>AR1533</u>	Graphite crucible for N/O, box of 1000 or bags of 100
905.200.020.001	AR20020	Graphite crucible for H2, box of 1000 or bags of 100
905.200.320.001	<u>AR936</u>	Spoon for weighing tungsten
905.200.380.001AN	<u>AR3818</u>	Ceramic crucibles 1000 pcs/box
905.200.430.001		O-ring P11 (FPM) reagent tube holder for upper and lower
905.200.550.001	<u>AR171</u>	Dehumidifier Anhydron 24-48 mesh .5 kg
905.200.560.001	<u>AR117</u>	Combustion tube transparent quartz glass 35mm x 148mm length
905.200.640.001		O-ring P4 (FPM) reagent tube holder filter holder
905.200.890.001		Reagent tube quartz 15mm dia. x 225mm length, purifier catalyst
905.200.920.001		Net filter 100 mesh SUS dia. 7mm for reagent tube
905.200.970.001	<u>AR2174</u>	Ascarite 20-30 mesh .5kg
905.200.970.001AN	<u>AR2174</u>	Ascarite Alphasolve 20-30 mesh .5kg
905.201.140.001	<u>AR8053</u>	Combustion boat ceramic No. 6a 500 pcs/box
905.201.190.001		Boat cover for No 6a boat 100 pcs/box
905.201.210.001	<u>AR022B</u>	Brush 34mm dia. cylindrical brush for dust removal
905.201.240.001		Flororubber packing for HF dust filter
905.201.250.001AN	<u>AR076</u>	Accelerator Particle-type tin .9kg

905.201.280.001		O-ring P35 (SI) heat-proof HF combustion tube seal
905.201.300.001		Work Coil for heating samples HF furnace
905.201.320.001		Oscillator tube 6T58RA HF furnace
905.201.340.001	<u>AR076</u>	Accelerator particle-type Tin 1kg
905.201.450.001		Reagent tube Pyrex 20m dia. x 150mm length
905.201.660.001	<u>AR9721</u>	Dust filter 1-6 micron quartz wool
905.201.790.001	<u>AR79001</u>	Graphite crucible, box of 1000 or bags of 200
905.201.810.001		Brushes for automatic cleaner 2 pcs
905.201.820.001		Crucible stand ceramic H63mm
905.201.830.001		Cylindrical filter stainless steel 2 micron furnace dust filter
905.201.850.001		Silicone tube for dust valve
905.201.890.001		Silicone tube 120mm for dust valve + fixing
905.201.910.001		Exhaust tube Ultrasoft flew for dust box
905.201.920.001		Mask filter 89-018-02 (5 pcs) white filter for dust box
905.201.930.001		Cylindrical filter stainless steel 2 micron furnace dust filter
905.201.940.001		Lower brush .3mm dia. for automatic cleaner
905.201.950.001		Upper brush .2mm dia. for automatic cleaner (optional for 320v and 220v)
905.202.030.001		O-ring P50 (heat resistant silicone rubber) for sealing the filter cap
905.202.040.001		Lower brush .2mm dia. for automatic cleaner (optional for 320v and 220v)
905.202.050.001		O-ring P8 (heat resistant silicone rubber) for sealing the reagent tube for purifier
905.202.060.001		O-ring JISW1516 No.29 (silicone rubber) semi transparent for sealing the furnace
905.202.070.001		O-ring P35 (FKM) for the combustion tube
905.202.080.001		X-ring X115A (FKM) for sealing the auto cleaner rods
905.202.090.001		Dust catcher, flame retarded fiber (black) filter for dust box
905.202.100.001	<u>AR5203HS</u>	Pedestal, ceramic ~24.5mm W x 50mm L, cup tip (Horiba style)
905.202.100.001	<u>AR5202LE</u>	Pedestal, ceramic ~ 24.5mm W x 50mm L, easy clean (LECO style)
905.202.110.001		Cylindrical filter stainless steel 2 micron furnace cylindrical filter

905.202.120.001		O-ring G35 (silicone rubber) cylindrical filter to seal both ends
905.202.140.001		Dust packing (FKM) for sealing lower part of furnace
905.202.140.001AN	<u>AR262</u>	Accelerator copper 20-60 mesh 1.4kg
905.202.150.001	<u>AR263</u>	Accelerator copper 20-60 mesh .5kg
905.202.200.001	<u>AR3818B</u>	Ceramic crucibles 500 pcs/box
905.2023.130.001		HEP filter (white) for dust box
905.900.330.001		Tygon tube for dust valve
905200.550.001AN	<u>AR171</u>	Dehumidifier Anhydron 24-48 mesh 454g
F002.720.200		Silicone tube for dust valve 1 meter
F0020.031.900		O-ring P10 (FPM) autocleaner for rod seal
F0020.433.100		O-ring G130 (FPM) dust box
F0020.945.500		O-ring P11 (SI) milk white HF furnace lower part flange
F0021.2942.800		O-ring S36 HF dust filter for both ends
G 012 09 20	<u>AR01046</u>	Platinum catalyst for EMIA purifier 40g
G881.14.315		Filter on side for ambient air 2 pcs
G881.62.54		Filter on side for ambient air
H0370.434.001		Exhaust tube Ultrasoft flew for dust box