

DUMATHERM #	ALPHA #	DESCRIPTION
14-0007	AR140007	Reduced Copper Stick 475g
14-0011	AR140011	HT Oxidation Catalyst 100g
14-0012	AR9721	Quartz Wool 28g
14-0013	AR171D	Magnesium Pechlorate 10g
14-0014	AR140014	LT Oxidation Catalyst 50g
14-0015	AR5000	Ceramic Ash Insert Package of 10
14-0016	ATD5216	Pressed Tin Capsule
14-0017	AR2186 or AR2286-200	Tin Foil Cups, -100 or 200 foils
14-0018	AR140018	Quartz Reaction Tube
14-0019	AR140019	Seals, bottom connection
14-0020	AR140020	Seals, top of reaction tubes
14-0026	AR140026	PrePacked Combustion Tube (filled with LT AND HT)
14-0027	AEB1002	Copper Oxide Wires
14-0029	AR5000	Ceramic Ash Insert Package of 10
14-0072	AR1618	Tweezer
14-0078	AR140078	O-ring for Autosampler (3)