


Product description:

PW-WIN software is designed for collecting measurement data from RADWAG weighing instruments and its statistics process.



Designed for scales:

- with standard communication frame

Functions:

- readout from instrument, record from measurements and saving data to a file
- possibility of manual and automatic registration of measurements
- registration and storing of measurement data in a form of file
- visualization of measurements data in a form of graph
- static processing of measurement data
- possibility of creating data from chosen measuring session
- printout of measuring data, graphs and statistics
- mainly used in laboratories

Computer requirements:

- processor 1,2 GHz
- free 500 MB on hard drive
- RAM 256 MB (recommended 512 MB)
- operation system Windows 2000/XP

Language versions:

Polish, German, English, French.

* NOTE! RADWAG does not sell software "in box". Boxes presented on our Web-site are of virtual characteristic only, and they are not accessible with purchased software. The customer receives purchased software on a CD in special protective cover and the user manual as a brochure.