

# PROFESSIONAL WEIGHING TABLE



Professional Weighing Table is intended for the use with microbalances and analytical balances. It incorporates an independent anti-vibration table with a stone table top. This solution eliminates possible vibrations and significantly shortens balance stabilization time.

Two lockable cabinets and pull-out keyboard tray facilitate installing PC computer with a keyboard, mouse, and UPS power supply. The raised monitor platform at the rear part of the table is designed for the installation of LCD monitor. The table features high functionality and ergonomics.

#### Unit contents

- weighing table
- anti-vibration table

