

# CEREAL DENSITY KIT



This kit is designed for determination of density of cereals. Measurement of cereal density in loose state is performed with application of precise balance series WPX 4500/G and density kit for cereals. Density determination result is calculated automatically by balance software (in accordance to tables containing cereal density).

It is possible to determine the density of following cereal: wheat, oats, barley, rye.

Manufactured of brass, the density kit is additional equipment to a balance. Its volume is 1 liter.

The cereal density kit allows for determination of density according to norm: PN-73/R-74007.

## CEREAL DENSITY KIT

- 1-base for measuring container
- 2-measuring container of 1 l volume
- 3-knife
- 4-ring
- 5-filling device
- 6-charging vessel

