

SUITCASE FOR SCALE



Special design providing protection for weighing instrument while transporting to location. Durable for shocks (aluminum construction with internal filling from polyurethane foam).

Manufactured in two versions:

- W1 for balances: WPX, PS xx/C2, PS xx/C1, WLL xx/A2, WLC xx/A1, WLC xx/A2 - with pan size: 195×195 mm.

- W2 for balances: WLC xx/C1/K, WTC xx/C1/K, WLL xx/C1/K - with pan size 290×360 mm.

