



PLEASE CONTACT US:
Tel: (416) 742-9068
Toll free: 1 (888) 742-9068
Fax: (416) 742-6817
www.romquest.com
sales@romquest.com

Code	Description	Pack
------	-------------	------

REFERENCE STANDARD MATERIAL

		C%	H%	N%	S%	O%
E11001	ACETANILIDE (1g)	71.09	6.71	10.36	-	11.84
E11002	4-FLUOROBENZOIC ACID (1g)	60.01	3.59	-	-	22.84
E11003	UREA (1g)	20.00	6.71	46.63	-	26.64
E11004	ANTHRACENE (1g)	94.34	5.66	-	-	-
E11005	EDTA (1g)	41.09	5.52	9.59	-	43.79
E11007	ATROPINE SULPHATE (1g)	58.77	7.25	4.03	4.62	25.33
E11009	L-CYSTINE (1g)	29.99	5.03	11.66	26.69	26.63
E11012	B-BOT (1g)	72.53	6.09	6.51	7.44	7.43
E11016	BENZOIC ACID (1g)	68.85	4.95	-	-	26.20
E11018	NICOTINIC ACID (1g)	58.54	4.09	11.38	-	-
E11024	1-Cl-2,4-DINITROBENZENE (1g)	35.58	1.49	13.83	-	31.60
E11028	SULPHANILAMIDE (1g)	41.85	4.68	16.26	18.62	18.58
E11052-B	LUBRICANT OIL (2g)	82.07	13.62	1.09	3.22	-

REFERENCE STANDARD MATERIALS

Accredited vs. NIST Std.

E11001-A	ACETANILIDE (1g)	71.09	6.71	10.36	-	11.84
E11003-A	UREA (1g)	-	-	46.63	-	-
E11009-A	L-CYSTINE (1g)	29.99	5.03	11.66	26.69	26.63

CERTIFIED REFERENCE STANDARD MATERIALS

E11030-C	B USH, BRANCHES & LEAVES (NS) (2g)	-	-	1.50	0.73	-
----------	------------------------------------	---	---	------	------	---

REFERENCE STANDARDS MATERIAL for SOILS

E11035-A	SYNTHETIC MIX. for SOIL #1 (5g)	3.500	-	0.216	-	-
E11036-A	SYNTHETIC MIX. for SOIL #2 (5g)	0.732	-	0.064	-	-
E11037-A	SYNTHETIC MIX. for SOIL #3 (5g)	4.401	-	0.262	-	-
E11038-A	SYNTHETIC MIX. for SOIL #4 (5g)	2.417	-	0.048	-	-
E11039-A	SYNTHETIC MIX. for SOIL #5 (5g)	0.141	-	0.021	-	-

Synthetic Mixtures #1 and #4 include contents of Organic and Inorganic Carbon differentiation:

		Total Organic Carbon	Total Inorganic Carbon
E11035-A	SYNTH. MIX. #1	2.613%	0.887%
E11038-A	SYNTH. MIX. #4	0.328%	2.089%