


Product description:

The software determines the percent of content of determined fractions (e.g. dimension of grains) in a sample.



Designed for scales:

- terminal TMC

Functions:

- defining the fractions
- data base on users and sieves
- definitions of analyses
- weighing of each fraction
- percent analysis of content of each fraction
- generating reports from measurements and from their analysis
- readout of measurements and their record
- record of measurement data in a form of a file from balance of measurements
- registration of control of measurements

Computer requirements:

- processor 1,2 GHz
- free 500 MB on hard drive
- RAM 256 MB (reccommended 512 MB)
- operation system Windows 2000/XP

* NOTE! RADWAG does not sell software "in box". Boxes presented on our Web-site are of virtual characteristic only, and they are not accessible with purchased software. The customer receives purchased software on a CD in special protective cover and the user manual as a brochure.